

Digital pH Meter Model : pHCal



Supplied with Combination Glass pH Electrode (0-70°C), Standard Buffer Tablets, Electrode Stand with Clamp, Dust Cover, Power Cord, Operational & Instruction Manual, Manufacturer Test Certificate & Warranty Card

Technical Specification :

- **pH Range** : 0-14 pH,
- **Resolution** : 0.01 pH
- **Temperature Range** : 0.0 to 100°C (Manual Compensation)
- **Display** : 3 1/2 Digit LED Display
- **Power Supply** : 230VAC \pm 10%, 50 Hz

Salient Features :

Calibration check facility & Calibration Error indication for 7.00 & 4.00 pH