

**Handheld pH/mV/Temperature Analyzer,  
Model :  $\mu$ pHCal<sub>Handheld</sub>**



Supplied with Combination Glass pH Electrode (0-70°C), Temperature Probe, Carry Bag, Standard Buffer Tablets, Two Bottle to fill Standard Buffer, Rechargeable Battery, DC Adaptor, Operational & Instruction Manual, Manufacturer Test Certificate & Warranty Card.

**Technical Specification :**

**pH Range :** -1.00 to 15.00 pH

**Resolution :** 0.01 pH

**Relative Accuracy :**  $\pm$  0.01 pH

**Backlit LCD Display :** 16 x 2 line Alphanumeric

**Auto Buffer Recognition :** 7.00 pH , 4.01/9.18 pH @25°C

**Temperature Range :** 0.0 to 130°C

**Power Supply :** Rechargeable Battery Operated supplied with DC Adaptor.

**Salient Features :** Auto Buffer Recognition & Simultaneous Display of pH, mV & Temperature, Auto / Manual Temp. Compensation.

**Error Message :**

- Bad Buffer / Electrode
- Temperature Probe
- pH over range