


pH/mV/°C Analyzer Model : μ pHCal₅ (Two Point Calibration)



Supplied with Combination Glass pH Electrode (0-70°C), Standard Buffer Tablets, Temperature Probe, Electrode Stand with Clamp, Dust Cover, DC Adaptor, Operational & Instruction Manual, Manufacturer Test Certificate & Warranty Card.

Technical Specification :

- pH Range : -1.00 to 15.00 pH
- Resolution : 0.01 pH
- Relative Accuracy : ± 0.01 pH
- mV : ± 2000.0
- Auto Buffer Recognition : 7.00, 4.01 / 9.18 @ 25°C
- Temperature Range : 0.0 to 130°C
- Backlit LCD Display : 16 x 2 line Alphanumeric
- Power Supply : DC Adaptor

Salient Features :

- Auto Buffer Recognition & Simultaneous Display of pH, mV & Temperature,
- Auto / Manual Temp. Compensation etc.

Error Message :

- Bad Buffer / Electrode ▪ Temperature Probe ▪ pH over range